INSPT0345/10.2025/V1

PT 0345 TEMBLAK kończyny górnej YOUNG

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INSTRUCTIONS FOR USE

NÁVOD NA POUŽITIE

INSTRUCIUNI DE UTILIZARE

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INSTRUCTIONS FOR USE

INTENDED USE:

Arm sling YOUNG is a reusable medical device intended to support treatment and rehabilitation according to the relevant indications.

INDICATIONS:

- post-traumatic (in situations that do not require full immobilisation or restriction of mobility)
- post-operative (in situations that do not require full immobilisation or restriction of mobility)
- limb relief, support for the upper extremity in a plaster dressing
- during rehabilitation

CONTRAINDICATIONS:

– no contraindications to using the product.

STORAGE:

Store the device in a dry place. Protect from high temperatures and direct sunlight.

TARGET GROUP:

The device is intended for people with the symptoms described in the recommendations.

PRECAUTIONS:

The device should be used only on intact skin. Do not apply directly to wounds, damaged or pathologically affected skin. The use of ointments, oils, creams and similar substances may adversely affect the properties of the device, reduce its durability and diminish its effectiveness. In the event of persistent or increased pain, allergic reaction, skin irritation or any other adverse effect related to the use of the device, discontinue use and consult a physician. The appropriate use of the device may require the advice of specialist medical staff. For hygiene reasons and to ensure effectiveness, the product is intended for use by a single patient only. Selection of the correct size will determine the correct functioning of the device. Each time the device is used, it should be checked for signs of wear or ageing. Only an undamaged device guarantees optimal performance. Unless otherwise recommended by your doctor, do not use the device continuously for 24 hours.

INFORMATION ON INITIAL PREPARATION OF THE PRODUCT:

It is recommended to wash the product before first use.

MAINTENANCE:

Fasten the Velcro before washing in order to avoid the device getting caught on something.

INSTRUCTIONS FOR USE:

Pass the tightening strap through the plastic frame, adjusting it from the inside depending on whether it is to be used on the left or right arm. Place your forearm in the arm sling pocket. Hook your thumb into the designated loop (optional) to prevent the device from moving up your arm. Pass the tightening strap through your healthy arm. Use the Velcro fastener to adjust the length of the tightening strap.

COMPOSITION PT 0345:

Polyester 66%, Polypropylene 27%, Polyamide 7%

PRODUCT ADJUSTMENT:

Size chart is available on the packaging.

DISPOSAL:

The product must be disposed of in accordance with applicable local waste management regulations.

Any serious incident related to the product must be reported to the manufacturer and competent authority of the Member State where the user or patient is resident.

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SIGVARIS S.A.